



## Cannabinoids Separation by Using mini LC-80 HPLC

### Instrumentation

HPLC condition:-

- HPLC : model mini LC-80 (Biometrics Technologies Inc.)
- Column : Raptor ARC-18, 2.7  $\mu\text{m}$ , length 150 mm and internal diameter 4.6 mm (Restek)
- Flowrate : 1.5 mL/min
- Detector : UV at 228 nm
- Isocratic system
- Mobile phase :  
A: 5 mM Ammonium formate in water, 0.1% formic acid  
B: Acetonitrile, 0.1% formic acid  
Mobile phase ratio A: B = 25 : 75
- Sample Injection volume: 20  $\mu\text{L}$



### Cannabinoids Standard

- Cannabinoids Standard (3 components) 1000  $\mu\text{g/mL}$ , Restek (Cannabinol (CBN), Cannabidiol (CBD),  $\Delta$ 9-Tetrahydrocannabinol ( $\Delta$ 9-THC))
- Cannabigerol (CBG) 1000  $\mu\text{g/mL}$ , Restek
- Tetrahydrocannabivarin (THCV) 1000  $\mu\text{g/mL}$ , Restek
- $\Delta$ 8-Tetrahydrocannabinol ( $\Delta$ 8-THC) 1000  $\mu\text{g/mL}$ , Restek

## Preparation of 100 ppm of cannabinoids solution (Stock solution)

Cannabinoids of 100 ppm was prepared by pipetting 100  $\mu\text{L}$  each standard bottom to the a vial and making it to the 1.0 mL with the mobile phase.

## Preparation of working cannabinoids standard solution

The working standard solution was prepared by pipetting 100, 200, 300, 400, 500  $\mu\text{L}$  stock solution into vials and then making it to the 1.0 mL with the mobile phase.

## Result and discussion




Figure 1. Chromatogram of 100 ppm cannabinoid standard 20  $\mu\text{L}$  Injection, 228 nm UV Detector

Figure 1 Shows the Chromatogram for the cannabinoids analysis. The chromatogram shows 6 peaks corresponding to in order of elution Cannabigerol (CBG) at 3.284 min, Cannabidiol (CBD) at 3.460 min, Tetrahydrocannabivarin (THCv) at 3.835 min, Cannabinol (CBN) at 5.33 min,  $\Delta^9$ -Tetrahydrocannabinol ( $\Delta^9$ -THC) at 6.817 min and  $\Delta^8$ -Tetrahydrocannabinol ( $\Delta^8$ -THC) at 7.042 min.




Figure 2. The calibration curve of mix cannabinoids standard solution A- Canabigerol (CBG), B- Cannabidiol (CBD) , C- Tetrahydrocannabivarin (THCV), D- Canabinol (CBN), E-  $\Delta$ 9-Tetrahydrocannabinol ( $\Delta$ 9-THC), F-  $\Delta$ 8-Tetrahydrocannabinol ( $\Delta$ 8-THC)

Figure 2 shows the calibration curve of cannabinoids which are prepared by using mix standard solution in concentration range of 10–50 ppm. The corresponding equation and correlation coefficient are shown.

## Conclusion

The compact mini LC-80 has been designed to give the chemist to have an easy-to-use HPLC for a wide range applications. The cannabinoids separation have been reported in this note showing off the application carried out in analytical laboratories can be easily and rapidly using mini HPLC.

